



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human IL-2 Monoclonal Antibody [5A8] Datasheet (Tracer, HRP conjugated)

Catalogue No	MO-C40032T
Product name	Mouse anti-human IL-2 monoclonal antibody [5A8]
Clone No	5A8
Product type	Primary antibody, Tracer
Description	Mouse monoclonal antibody to human IL-2, Horseradish Peroxidase (HRP) conjugated.
Target protein	Human IL-2
Immunogen	Recombinant Human IL-2
Reactivity	Reactive to the recombinant Human IL-2
Cross-reactivity	Cross-reactivity was not identified
Clonality	Monoclonal
Source	Mouse
Myeloma	Sp2/0-Ag14
Subclass	IgG2b
Light Chain	Kappa
Formulation	Lyophilized in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended to adjust the final concentration to 1.00 mg/ml
Applications	ELISA: HRP conjugated Clone 5A8 can be used as detection antibody to detect recombinant human IL-2 in ELISA when clone 9F9 is used as capture antibody.
Research	Cytokine
Storage	Store at -20° C. Avoid repeated freeze and thaw cycles
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)