



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Human Umbilical Cord Serum Datasheet

Catalogue No	S30010
Product name	Human umbilical cord serum
Description	Human umbilical cord serum is prepared from new born placenta umbilical cord which is non-reactive to hepatitis B surface antigen, HIV-1 antibodies. The AFP content was \geq 0.05mg/ml. Please note: For cell culture and other applications that require sterile condition, sterile filtration is required.
Source	Human blood
Storage	Short Term: 4°C for a few days. Long Term: -20°C.
Format	Liquid
Formulation	Liquid
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(02)