



## **ANOGEN - A Division of YES Biotech Laboratories Ltd.**

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

### **Mouse anti-Hepatitis B Surface Antigen (HBsAg) Monoclonal Antibody [A6] Datasheet**

<b>Catalog No</b>	MO-I40013B
<b>Product Name</b>	Mouse anti-Hepatitis B Surface Antigen (HBsAg) monoclonal antibody [A6]
<b>Clone number</b>	A6
<b>Product type</b>	Primary antibodies
<b>Description</b>	Mouse monoclonal antibody to Hepatitis B Surface Antigen (HBsAg), Protein G affinity purified
<b>Target protein</b>	HBsAg
<b>Reactivity</b>	Human, others not tested
<b>Clonality</b>	Monoclonal
<b>Raised in</b>	Mouse
<b>Subclass</b>	IgG1
<b>Light chain type</b>	Kappa
<b>Formulation</b>	Lyophilized from a solution in 0.01M PBS, pH 7.2
<b>Reconstitution Solution</b>	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
<b>Applications</b>	ELISA
<b>Research areas</b>	Virology
<b>Storage instructions</b>	Store at -20 <sup>0</sup> C. After reconstitution, aliquot and store at -20 <sup>0</sup> C.
<b>References</b>	<i>If you publish research using this product, please let us know so that we can cite the reference on this datasheet. In return we will send you a reward of your choice.</i>

**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**

S7.5(01)