



## **ANOGEN - A Division of YES Biotech Laboratories Ltd.**

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

### **Rat anti-Mouse CD8 Monoclonal Antibody [JXYT8] Datasheet**

<b>Catalogue No</b>	MO-L400027
<b>Product name</b>	Rat anti-mouse CD8 monoclonal antibody [JXYT8]
<b>Clone No</b>	JXYT8
<b>Product type</b>	Primary antibodies
<b>Description</b>	Rat monoclonal antibody (MAb) to mouse CD8
<b>Target protein</b>	Murine T-cell surface antigen CD8 marker
<b>Immunogen</b>	Mouse cytolytic T-lymphocytes
<b>Specificity</b>	Reactive with mouse T-cell surface antigen CD8 marker
<b>Reactivity</b>	Mouse, others not tested.
<b>Clonality</b>	Monoclonal
<b>Source</b>	Rat
<b>Subclass</b>	IgM
<b>Formulation</b>	Lyophilized from a solution in 0.01M PBS pH7.2
<b>Reconstitution</b>	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
<b>Storage</b>	Store at -20° C
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**