



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti- TRICHOSANTHIN (TCST) Monoclonal Antibody Datasheet [TCST1]

Catalogue No	MO-M40011A
Product name	Mouse TRICHOSANTHIN (TCST) monoclonal antibody [TCST1]
Clone No	TCST1
Product type	Primary antibodies
Description	Mouse monoclonal antibody (MAb) to TCST Protein G affinity purified
Target protein	Human TRICHOSANTHIN (TCST)
Immunogen	Highly purified TCST
Specificity	Reactive with TCST
Reactivity	TCST , others not tested
Cross-reactivity	unknown
Clonality	Monoclonal
Source	Mouse
Subclass	IgG1
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
Applications	ELISA
Research	Family plan
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)