



**ANOGEN** - A Division of YES Biotech Laboratories Ltd.  
2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

**Mouse anti-Human STAT5b Monoclonal Antibody [A2H6] Datasheet**

<b>Catalogue No</b>	MO-P40001B
<b>Product name</b>	Mouse anti-Human STAT5b (p-T694) monoclonal antibody.
<b>Clone No</b>	A2H6
<b>Product type</b>	Primary antibody.
<b>Description</b>	Mouse monoclonal antibody to human STAT5b, Protein G affinity purified.
<b>Target protein</b>	Human STAT5b
<b>Immunogen</b>	Peptide of a phosphorylated fragment of STAT5b (p-T694), conjugated with KLH
<b>Specificity</b>	Reactive to the immunogen and peptide-BSA conjugating.
<b>Reactivity</b>	Human, others not tested.
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mouse
<b>Subclass</b>	IgG2b
<b>Light Chain</b>	Kappa
<b>Formulation</b>	Lyophilized form from a solution of 0.01M PBS, pH7.0
<b>Reconstitution</b>	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
<b>Applications</b>	Indirect ELISA. The suggested antigen coating quantity is 1 ug/ml of peptide/or peptide conjugating in pH 9.5 carbonate buffer, 100 ul/well .
<b>Research Area</b>	Phosphorylation, metabolism.
<b>Storage</b>	Store at -20 <sup>o</sup> C. After reconstitution, aliquot and store at -20 <sup>o</sup> C; avoid thaw and freeze cyclically.
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**