




**Mouse anti-GST Tag Monoclonal Antibody [GST.B6] Datasheet**

<b>Catalogue No</b>	MO-H40003A
<b>Product name</b>	Mouse anti-GST Tag Monoclonal Antibody [GST.B6]
<b>Clone No</b>	GST.B6
<b>Product type</b>	Purified primary antibodies
<b>Description</b>	Mouse monoclonal antibody against GST (glutathione S-transferase) tag.
<b>Target protein</b>	GST and GST fusion protein
<b>Immunogen</b>	A synthetic peptide from GST coupled to KLH.
<b>Specificity</b>	This monoclonal antibody binds to GST-tagged recombinant protein.
<b>Cross-reactivity</b>	Not identified
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mouse
<b>Subclass</b>	IgG2a
<b>Light Chain</b>	Kappa
<b>Formulation</b>	Liquid, in pH7.4 PBS with 0.05% sodium azide.
<b>Storage</b>	Aliquot and store at -20°C for long term storage. Avoid repeated freeze and thaw cycles.
<b>Applications</b>	Western blotting (1:1,000). Immunoprecipitation and Immunostaining (1:200)
<b>Research</b>	Life science
<b>Storage</b>	Store at -20° C
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

**Experimental data**



**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**

S7.5(01)