




ANOGEN - A Division of YES Biotech Laboratories Ltd.
2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023


Mouse anti-His Tag Monoclonal Antibody [HIS.H8] Datasheet

Catalogue No	MO-H40005A
Product name	Mouse anti-His Tag Monoclonal Antibody [HIS.H8]
Clone No	HIS.H8
Product type	Purified primary antibodies
Description	Mouse monoclonal antibody against synthetic 6 x Histidine tag.
Target protein	His-tagged recombinant protein
Immunogen	Synthetic 6 x Histidine peptide
Specificity	This monoclonal antibody binds to His-tagged recombinant protein.
Cross-reactivity	Not identified
Clonality	Monoclonal
Source	Mouse
Subclass	IgG2b
Light Chain	Kappa
Formulation	Liquid, in pH 7.2 0.01M PBS with 0.05% sodium azide.
Storage	Aliquot and store at -20°C for long term storage. Avoid repeated freeze and thaw cycles.
Applications	Western blotting (1:1,000). Immunoprecipitation and Immunostaining (1:2,000)
Research	Life science
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

Experimental data



His Western Blot analysis of 293 cells transfected with His-tagged vector (2) and untransfected control (1).



Immunofluorescence staining of His-tag fusion protein in 293 cells. A: untransfected control, B: transfected

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)