




ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-His Tag Monoclonal Antibody [10A12] Datasheet

Catalogue No	MO-H40005B
Product name	Mouse anti-His Tag Monoclonal Antibody [10A12]
Clone No	10A12
Product type	Purified primary antibodies
Description	Mouse monoclonal antibody against synthetic 6 x Histidine tag
Target protein	His-tagged recombinant protein
Immunogen	Synthetic 6 x Histidine peptide conjugated with KLH
Specificity	This monoclonal antibody binds to His-tagged recombinant protein, but not with KLH
Cross-reactivity	Not identified
Clonality	Monoclonal
Source	Mouse
Subclass	IgG1
Light Chain	Kappa
Formulation	Lyophilized in pH 7.2 0.01M PBS.
Reconstitution	Reconstitute with appropriate volume of dH ₂ O to reach a final concentration of 1mg/ml
Applications	Western blotting
Research	Life science
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

Experimental data: Western blotting using anti-His antibodies (HisH8 and 10A12) to detect His-tagged HIV p24



This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)